




Enfermedades Infecciosas y Microbiología Clínica

www.elsevier.es/eimc



Brief report

Comparison of three procedures for the rapid identification of bacteraemia-causing microorganisms. Evaluation of their effectiveness and applicability to microbiology laboratories[☆]

Oriol Martín-Pujol, Tomas Tosco-Nuñez*, Isabel de Miguel-Martinez

Servicio de Microbiología, Hospital Universitario Insular de Gran Canaria, Las Palmas de Gran Canaria, Spain

ARTICLE INFO

Article history:

Received 5 March 2018

Accepted 1 June 2018

Available online 2 April 2019

Keywords:

Blood cultures

Bacteraemia

Identification

Rapid methods

MALDI-TOF

ABSTRACT

Introduction: Three procedures for rapid identification of microorganisms in positive blood cultures were evaluated.

Methods: We performed two methods based on direct extraction from a blood culture: Sepsityper[®] (Bruker Daltonics) (ST) and a non-commercial saponin method (MCS), and another method consisting of a short incubation subculture (SIC). Identification values obtained by spectrometry Matrix-Assisted Laser Desorption Ionization-Time of Flight (EM MALDI-TOF) were compared by applying the manufacturer's interpretation criteria and corrected cut-off points.

Results: According to the manufacturer, 65.8%, 45.8% and 57.4% of microorganisms were identified at the species level by using ST, MCS and SIC, respectively. When applying corrected cut-off points, the values increased to 92.3%, 80.6% and 85.2%, respectively. ST offered significantly better results than MCS, and no significant differences were found between ST and SIC, except for with respect to yeast.

Conclusions: Better identification rates were obtained by using ST and SIC, which are easily applicable in any laboratory.

© 2018 Elsevier España, S.L.U. and Sociedad Española de Enfermedades Infecciosas y Microbiología Clínica. All rights reserved.

Comparación de tres procedimientos para la identificación rápida de microorganismos causantes de bacteriemias. Evaluación de su eficacia y aplicabilidad en el laboratorio de microbiología

RESUMEN

Introducción: Evaluamos tres procedimientos de identificación rápida de microorganismos a partir de hemocultivos positivos.

Métodos: Aplicamos dos métodos basados en la extracción directa desde el frasco de hemocultivo: Sepsityper[®] (Bruker Daltonics) (ST) y un método casero con saponina (MCS), y un tercer método basado en un subcultivo con incubación corta (SIC). Se comparan las identificaciones por espectrometría de masas Matrix-Assisted Laser Desorption Ionization-Time of Flight (EM-MALDI-TOF) aplicando los criterios de interpretación del fabricante y los puntos de corte corregidos (PCC).

Palabras clave:

Hemocultivos

Bacteriemia

Identificación

Métodos rápidos

MALDI-TOF

DOI of original article: <https://doi.org/10.1016/j.eimc.2018.06.018>

[☆] Please cite this article as: Martín-Pujol O, Tosco-Nuñez T, de Miguel-Martinez I. Comparación de tres procedimientos para la identificación rápida de microorganismos causantes de bacteriemias. Evaluación de su eficacia y aplicabilidad en el laboratorio de microbiología. *Enferm Infecc Microbiol Clin.* 2019;37:319–323.

* Corresponding author.

E-mail address: tomastn@hotmail.com (T. Tosco-Nuñez).

Resultados: Aplicando los criterios del fabricante se identificaron a nivel de especie el 65,8%, el 45,8% y el 57,4% con ST, MCS y SIC, respectivamente. Aplicando los PCC, estos resultados fueron del 92,3%, del 80,6%, y del 85,2%, respectivamente. La identificación con ST fue significativamente mejor que el MCS. ST y SIC no mostraron diferencias significativas, excepto en levaduras.

Conclusiones: ST y SIC obtienen buenas tasas de identificación y pueden integrarse fácilmente en cualquier laboratorio.

© 2018 Elsevier España, S.L.U. y Sociedad Española de Enfermedades Infecciosas y Microbiología Clínica. Todos los derechos reservados.

Introduction

The aim of an effective empirical antibiotic therapy is to improve patient prognosis while reducing length of hospital stay and associated healthcare costs.¹

Köck et al.² found that clinicians changed the empirical treatment in 20% of their patients after identifying the species, versus 8% when only microscopy results were available. In light of these and similar data, new mass spectrometry identification methods like Matrix-Assisted Laser Desorption/Ionisation-Time of Flight (MALDI-TOF) have been developed based on positive blood culture bottles, both direct methods using commercial or in-house extraction techniques,^{3–10} as well as indirect methods using short-term subcultures.^{11–14} To interpret the scores obtained by MALDI-TOF, Kohlmann et al.¹¹ published some corrected cut-off points referring to Gram-positive and non-Gram-positive microorganisms.

The primary objective of this study was to evaluate the efficacy of three procedures for the rapid identification of microorganisms in positive blood cultures and their applicability in routine laboratory practice. The second objective was to compare the MALDI-TOF scores by applying the manufacturer's interpretation criteria as well as the corrected cut-off points proposed by Kohlmann et al. to assess their impact on the final identification of the microorganisms.

Material and methods

This comparative, observational study conducted between 2014 and 2016 evaluated three methods for the rapid identification of pathogens in positive blood cultures.

One bottle per patient and episode of bacteraemia was considered. The BacT/ALERT system (bioMérieux, Marcy l'Etoile, France) was used. The positive bottles were processed following standard protocols. After 24 h of incubation, the isolated microorganisms were identified by MALDI-TOF (Bruker Daltonics, Bremen, Germany) and by automated biochemical tests (MicroScan WalkAway, Beckman Coulter, USA). The identity of *Streptococcus pneumoniae* and *viridans* streptococci was confirmed by sensitivity to optochin and solubility in bile. Positive blood cultures with Gram stains that revealed polymicrobial aetiology were excluded.

Three aliquots were obtained from the first positive bottle of each patient studied, which were used to evaluate the following procedures:

- 1) Sepsityper[®] Kit (ST) (Bruker Daltonics, Bremen, Germany). 1 ml of blood culture was processed following the manufacturer's specifications and extraction was performed using ethanol and formic acid. 1 µl of supernatant was then placed on a metal plate.
- 2) In-house saponin method. Based on the protocol modification by Martiny et al.⁵ 1 ml was extracted from the blood culture bottle and 500 µl of a 5% saponin solution were added. After stirring in a vortex for 10 s, it was left to rest for 5 min at room temperature. It was washed twice with Milli-Q water. Ethanol-formic acid extraction was then performed and 1 µl of the supernatant was placed on a metal plate.

- 3) Short-term subculture (STS). A chocolate agar medium (Becton Dickinson, Franklin Lakes, NJ, USA) was spiked with 10 drops (≈500 µl) of blood culture and incubated in an oven at 37 °C, 5% CO₂ for 3.5 h. After this time had elapsed, a sterile inoculation loop was dragged several times across the surface of the culture and placed on a metal plate. Direct extraction on the plate with formic acid ≥96% was then performed.

Identification by MALDI-TOF

1 µl of stock solution (α-cyano-4-hydroxycinnamic acid) was added to the inoculum on the plate and left to dry at room temperature. The plate was then inserted in the MALDI-TOF system. The Bacterial Test Standard was used to calibrate the equipment, while the Microflex LT mass spectrometer and the software Biotyper 3.1 (Bruker Daltonics, Bremen, Germany) were used to analyse each sample.

The manufacturer's interpretation criteria defined values ≥2.3 as having a high identification probability, 2–2.2 as probable species identification and 1.7–1.9 as probable genus identification, while values <1.7 were interpreted as unreliable. According to the corrected cut-off points published by Kohlmann et al.,¹¹ correct genus and species identification was assumed with a score ≥1.5 for Gram-positive microorganisms and ≥1.7 for non-Gram-positive microorganisms. As an additional criterion, the first three identification options had to be identical or, failing that, the first two had to be identical with a difference greater than 0.3 with the third. Scores inconsistent with the above were considered to be unreliable. The direct identifications by MALDI-TOF that did not match the definitive identifications of the traditional subculture were considered to be incorrect.

Statistical analysis

The Student's *t* test with binomial distribution was used to compare the methods used. The total number of correct and unreliable/incorrect identifications were analysed for each method and an analysis by subgroup (Gram-positive microorganisms, Gram-negative microorganisms and yeasts) was also performed.

Results

In total, 155 blood culture bottles were analysed (90 aerobic and 65 anaerobic): 89 (57.4%) isolated from Gram-positive microorganisms, 53 (34.2%) from Gram-negative microorganisms and 13 (8.4%) from yeasts. Thirty-three (33) different species were identified. The overall results by microorganism and method, applying the corrected cut-off points,¹¹ are shown in Table 1. Table 2 presents the results of each method by interpretation criterion.

With the exception of two strains that were initially identified as *viridans* group streptococci before being ultimately identified as *S. pneumoniae* by manual testing, all the identifications obtained by the three study methods matched the identification of the traditional subculture.

Table 1
MALDI-TOF identification results by microorganism and methods used. The results were interpreted by applying the Kohlmann criteria.

ID after traditional subculture	Total N	Sepsityper®		In-house saponin method		Short-term subculture	
		Correct ID (species)	UID/No ID	Correct ID (species)	UID/No ID	Correct ID (species)	UID/No ID
Gram-negative bacteria	53	51 (96.2%)	2 (3.8%)	46 (86.8%)	7 (13.2%)	49 (92.5%)	4 (7.5%)
<i>Enterobacteria</i>	38	38 (100%)		36 (94.7%)	2 (5.3%)	36 (94.7%)	2 (5.3%)
<i>Citrobacter freundii</i>	1	1		1		1	
<i>Escherichia coli</i>	21	21		20	1	19	2
<i>Klebsiella pneumoniae</i>	6	6		5	1	6	
<i>Klebsiella oxytoca</i>	1	1		1		1	
<i>Morganella morganii</i>	2	2		2		2	
<i>Proteus mirabilis</i>	4	4		4		4	
<i>Salmonella</i> sp.	1	1		1		1	
<i>Serratia marcescens</i>	2	2		2		2	
<i>Non-fermenters</i>	13	13 (100%)		10 (76.9%)	3 (23.1%)	12 (92.3%)	1 (7.7%)
<i>Pseudomonas aeruginosa</i>	12	12		10		12	
<i>Pseudomonas oryzihabitans</i>	1	1			1		1
<i>Other</i>	2		2		2	1	1
<i>Haemophilus influenzae</i>	1		1		1	1	
<i>Helicobacter cinaedi</i>	1		1		1		1
Gram-positive bacteria	89	84 (94.4%)	5 (5.6%)	77 (86.5%)	12 (13.5%)	82 (92.1%)	7 (7.9%)
<i>Gram-positive cocci</i>	87	83 (95.4%)	4 (4.6%)	77 (88.5%)	10 (11.5%)	82 (94.3%)	5 (5.7%)
<i>Staphylococcus</i> spp.	56	56 (100%)		55 (98.2%)	1 (1.8%)	55 (98.2%)	1 (1.8%)
<i>Staphylococcus aureus</i>	16	16		16		16	
Coagulase-negative	40	40		39	1	39	1
<i>Staphylococcus capitis</i>	3	2		2		2	
<i>Staphylococcus epidermidis</i>	23	23		23		23	
<i>Staphylococcus haemolyticus</i>	1	1		1		1	
<i>Staphylococcus hominis</i>	12	12		11	1	11	1
<i>Staphylococcus lugdunensis</i>	1	1		1		1	
<i>Enterococcus</i> spp.	16	16 (100%)		12 (75%)	4 (25%)	16 (100%)	
<i>Enterococcus faecalis</i>	8	8		8		8	
<i>Enterococcus faecium</i>	8	8		4	4	8	
<i>Streptococcus</i> spp.	15	11 (73.3%)	4 (26.7%)	10 (66.7%)	5 (33.3%)	11 (73.3%)	4 (26.7%)
<i>Streptococcus agalactiae</i>	2	2		2		2	
<i>Streptococcus anginosus</i>	1	1			1	1	
<i>Streptococcus canis</i>	1	1		1		1	
<i>Streptococcus gallolyticus</i>	2	2		1	1	2	
<i>Streptococcus mitis</i>	1	1		1		1	
<i>Streptococcus oralis</i>	2	2		1	1	2	
<i>Streptococcus pneumoniae</i>	6	2	4	4	2	2	4
<i>Gram-positive bacilli</i>	2	1	1		2		2
<i>Lactobacillus rhamnosus</i>	1		1		1		1
<i>Propionibacterium acnes</i>	1	1			1		1
Yeasts	13	8 (61.5%)	5 (38.5%)	2 (15.4%)	11 (84.6%)	1 (7.7%)	12 (92.3%)
<i>Candida albicans</i>	3	1	2		3	1	2
<i>Candida glabrata</i>	1		1		1		1
<i>Candida parapsilosis</i>	8	6	2	2	6		8
<i>Candida tropicalis</i>	1	1			1		1
Total isolates	155	143 (92.3%)	12 (7.7%)	125 (80.6%)	30 (19.4%)	132 (85.2%)	23 (14.8%)

ID: identification; UID/No ID: unreliable identification/no identification.

The microorganism identification results from the ST were significantly better ($p < 0.05$) than the in-house saponin method. No significant differences were found for either Gram-negative ($p = 0.094$) or Gram-positive microorganisms ($p = 0.234$) in the ST and STS. However, the yeast analysis yielded a $p < 0.05$ in favour of ST. Finally, significant differences ($p < 0.05$) were found in the identification of both Gram-negative and Gram-positive microorganisms after comparing the in-house saponin method with the STS, with the latter producing better results.

Discussion

We analysed the effectiveness of three methods for the rapid identification of microorganisms in positive blood cultures. When the manufacturer's criteria were applied, the ST (65.8%) and the STS (57.4%) obtained the best species identification results. Significant differences were only found in the yeast identification parameter, with the commercial method yielding the best results. The results

of the in-house saponin method were significantly inferior (45.8%) to the other methods.

When the corrected cut-off points¹¹ were applied, the species identification percentages increased: ST (92.3%), STS (85.2%) and in-house saponin method (80.6%), which confirms that the manufacturer's cut-off points can be modified without losing specificity.^{12,15}

The results for the ST are similar to those already published, with correct species identification ranging from 59.5% to 82.8% when the manufacturer's criteria are applied.^{3–6} Our STS results are similar to those published in the literature.^{11,13,14} In terms of the in-house saponin method, the protocol amendments proposed by Martiny et al.⁵ yield better results when the corrected cut-off points¹¹ are applied.

When the corrected cut-off points¹¹ are applied, the Gram-negative bacteria species identification results of both the ST and the STS are excellent (96.2% and 92.5%, respectively). The overall Gram-negative bacteria identification results for the in-house saponin method were statistically inferior ($p < 0.05$), while the

Table 2Results obtained by Sepsityper[®], the in-house saponin method and the short-term subculture as per the interpretation criteria applied to the scores obtained by MALDI-TOF mass spectrometry.

Microorganisms isolated	N	Sepsityper [®]			In-house saponin method			Short-term subculture		
		Correct ID (genus) Manufacturer	Correct ID (species) Manufacturer	Correct ID (species) Corrected cut-offs	Correct ID (genus) Manufacturer	Correct ID (species) Manufacturer	Correct ID (species) Corrected cut-offs	Correct ID (genus) Manufacturer	Correct ID (species) Manufacturer	Correct ID (species) corrected cut-offs
<i>Gram-negative bacilli</i>	53	51 (96.2%)	47 (88.7%)	51 (96.2%)	47 (88.7%)	32 (60.4%)	46 (86.8%)	49 (92.4%)	43 (81.1%)	49 (92.5%)
Enterobacteria	38	28 (73.7%)	27 (71.1%)	38 (100%)	36 (94.7%)	27 (71.1%)	36 (94.7%)	36 (94.7%)	33 (86.8%)	36 (94.7%)
Non-fermenters	13	13 (100%)	10 (76.9%)	13 (100%)	11 (84.6%)	5 (38.5%)	10 (76.9%)	12 (92.3%)	9 (69.2%)	12 (92.3%)
Other	2	0	0	0	0	0	1	1	1	1
<i>Gram-positive cocci</i>	87	82 (94.2%)	54 (62.1%)	83 (95.4%)	64 (73.6%)	39 (44.8%)	77 (88.5%)	79 (90.8%)	46 (52.9%)	82 (92.1%)
Staphylococci	56	56 (100%)	34 (60.7%)	56 (100%)	46 (82.1%)	31 (55.4%)	55 (98.2%)	55 (98.2%)	27 (48.2%)	55 (98.2%)
<i>Staphylococcus aureus</i>	16	16	15	16	16	13	16	16	16	16
CNS	40	40	19	40	40	18	39	39	11	39
Enterococcus	16	15 (93.8%)	12 (75.0%)	16 (100%)	11 (68.7%)	8 (50.0%)	12 (75%)	15 (93.7%)	13 (81.3%)	16 (100%)
Streptococcus	15	11 (73.3%)	8 (53.3%)	11 (73.3%)	7 (46.7%)	0	10 (66.7%)	9 (60.0%)	6 (40.0%)	11 (73.3%)
<i>Gram-positive bacilli</i>	2	1 (50%)	0	1 (50%)	0	0	0	0	0	0
Yeasts	13	8 (61.5%)	1 (7.7%)	8 (61.5%)	2 (15.4%)	0	2 (15.4%)	1 (7.7%)	0	1 (7.7%)
Total isolates	155	142 (91.6%)	102 (65.8%)	143 (92.3%)	113 (72.9%)	71 (45.8%)	125 (80.6%)	129 (83.2%)	89 (57.4%)	132 (85.2%)

CNS: coagulase-negative staphylococci; ID: identification.

Enterobacteria results were comparable with the other two methods.

The ST and STS methods also obtained optimal results for Gram-positive bacteria when the corrected cut-off points¹¹ were applied (94.4% and 92.1%, respectively), with better identification of staphylococci and enterococci than streptococci. Not being able to distinguish between *S. pneumoniae* and *viridans* group streptococci is relatively common due to the similarity of their protein spectra. It is noteworthy that all strains of *Staphylococcus aureus* were correctly identified by the three methods and that no coagulase-negative staphylococcus was erroneously identified as *S. aureus*.

The yeast identification results were inferior to the rest, with the ST (61.5% applying the corrected cut-off points¹¹) statistically superior to the other methods.

In terms of the applicability of these techniques in laboratory practice, the ST is rapid (20 min) but entails a longer technical processing time and higher financial costs. The STS takes longer (3.5 h) but requires less technical processing and is cheaper, which means it could easily be integrated into routine laboratory practice.

Therefore, since the identification rates of the ST and STS methods are similar (except for yeasts), both techniques could be used in combination depending on the time to positivity of the blood culture and laboratory logistics, so as to maintain continuous clinical information with sustainable resource use.

Conflicts of interest

The authors declare that they have no conflicts of interest.

References

- Delpont JA, Strikwerda A, Armstrong A, Schaus D, John M. MALDI-ToF short incubation identification from blood cultures is associated with reduced length of hospitalization and a decrease in bacteremia associated mortality. *Eur J Clin Microbiol Infect Dis*. 2017;36:1181–6.
- Köck R, Wüllenweber J, Horn D, Lanckohr C, Becker K, Idelevich EA. Implementation of short incubation MALDI-TOF MS identification from positive blood cultures in routine diagnostics and effects on empiric antimicrobial therapy. *Antimicrob Resist Infect Control*. 2017;6:12.
- Kok J, Thomas LC, Olma T, Chen SCA, Iredell JR. Identification of bacteria in blood culture broths using matrix-assisted laser desorption-ionization Sepsityper™ and time of flight mass spectrometry. *PLoS ONE*. 2011;6:2–8.
- Loonen AJM, Jansz AR, Stalpers J, Wolffs PFG, van den Brule AJC. An evaluation of three processing methods and the effect of reduced culture times for faster direct identification of pathogens from BacT/ALERT blood cultures by MALDI-TOF MS. *Eur J Clin Microbiol Infect Dis*. 2012;31:1575–83.
- Martiny D, Dediste A, Vandenberg O, Vandenberg O. Comparison of an in-house method and the commercial Sepsityper™ kit for bacterial identification directly from positive blood culture broths by matrix-assisted laser desorption-ionisation time-of-flight mass spectrometry. *Eur J Clin Microbiol Infect Dis*. 2012;31:2269–81.
- Buchan BW, Riebe KM, Ledebner NA. Comparison of the MALDI biotyper system using Sepsityper specimen processing to routine microbiological methods for identification of bacteria from positive blood culture bottles. *J Clin Microbiol*. 2012;50:346–52.
- Meex C, Neuville F, Descy J, Huynen P, Hayette MP, de Mol P, et al. Direct identification of bacteria from BacT/ALERT anaerobic positive blood cultures by MALDI-TOF MS: MALDI Sepsityper kit versus an in-house saponin method for bacterial extraction. *J Med Microbiol*. 2012;61:1511–6.
- Yonetani S, Ohnishi H, Ohkusu K, Matsumoto T, Watanabe T. Direct identification of microorganisms from positive blood cultures by MALDI-TOF MS using an in-house saponin method. *Int J Infect Dis*. 2016;52:37–42.
- Tanner H, Evans JT, Gossain S, Hussain A. Evaluation of three sample preparation methods for the direct identification of bacteria in positive blood cultures by MALDI-TOF. *BMC Res Notes*. 2017;10:48.
- Chien J-Y, Lee T-F, Du S-H, Teng S-H, Liao C-H, Sheng W-H, et al. Applicability of an in-house saponin-based extraction method in Bruker Biotyper matrix-assisted laser desorption/ionization time-of-flight mass spectrometry system for identification of bacterial and fungal species in positively flagged blood cultures. *Front Microbiol*. 2016;7:1432.
- Kohlmann R, Hoffmann A, Geis G, Gatermann S. MALDI-TOF mass spectrometry following short incubation on a solid medium is a valuable tool for rapid pathogen identification from positive blood cultures. *Int J Med Microbiol*. 2015;305:469–79.
- Bhatti MM, Boonlayangoor S, Beavis KG, Tesic V. Rapid identification of positive blood cultures by matrix-assisted laser desorption ionization-time of flight mass spectrometry using prewarmed agar plates. *J Clin Microbiol*. 2014;52:4334–8.
- Idelevich EA, Schüle I, Grünastel B, Wüllenweber J, Peters G, Becker K. Rapid identification of microorganisms from positive blood cultures by MALDI-TOF mass spectrometry subsequent to very short-term incubation on solid medium. *Clin Microbiol Infect*. 2014;20:1001–6.
- Ballester-Téllez M, Recacha E, de Cueto M, Pascual A. Identificación y determinación de sensibilidad a antibióticos de aislados de hemocultivos a partir de subcultivos de corta incubación. *Enferm Infecc Microbiol Clin*. 2017;35:582–5.
- Hoyos-Mallecot Y, Miranda-Casas C, Cabrera-Alvargonzalez JJ, Gómez-Camarasa C, Pérez-Ramírez MD, Navarro-Marí JM. Identificación bacteriana directamente del hemocultivo mediante una técnica rápida de espectrometría de masas *Matrix-Assisted Laser Desorption-Ionisation Time-of-Flight*. *Enferm Infecc Microbiol Clin*. 2013;31:152–5.