



DIDACTIC

## Photoprocesses in chemical education. Key experiments for core concepts



Michael Walter Tausch\*, Nico Meuter, Sebastian Spinnen

*Department of Chemistry and Chemical Education, University of Wuppertal, Germany*

Received 20 January 2017; accepted 15 March 2017

Available online 17 April 2017

### KEYWORDS

Fluorescence;  
Phosphorescence;  
Ground states and excited states;  
Energy level model;  
Stokes shift;  
Down conversion;  
Photochromism;  
Solvatochromism;  
Aggregation Induced Emission (AIE)

**Abstract** Due to the “big five” challenges in the 21st Century related to the terms of energy, food, water, mobility and climate, mankind has to replace energy from fossil fuels step by step by renewable energies. Between them, solar energy is definitely the most abundant and clean. This is a major, but not the only reason for the inclusion of photoprocesses into chemical education. In this article we present an experimental approach to the basic concept suitable for teaching all phenomena involving light in a reasonable approximation. Updated experiments demonstrating the down conversion of light by fluorescence and room-temperature phosphorescence as well as aggregation induced emission and rigidification induced luminescence are described. In order to use these experiments for teaching basic concepts in photochemistry, a simple version of an energy level model, and an extended version, containing vibrational states, are proposed.

© 2017 Universidad Nacional Autónoma de México, Facultad de Química. This is an open access article under the CC BY-NC-ND license (<http://creativecommons.org/licenses/by-nc-nd/4.0/>).

### PALABRAS CLAVE

Fluorescencia;  
Fosforescencia;  
Estados base y estados de excitación;  
Modelo de nivel de energía;

**Fotoprocesos en la enseñanza de la química. Experimentos clave para conceptos básicos**

**Resumen** Debido a los «cinco grandes» retos del siglo XXI en relación con los términos de energía, alimentos, agua, movilidad y clima, la humanidad tiene que sustituir la energía de combustibles fósiles por las energías renovables paso a paso. Entre ellas, la energía solar es, sin ninguna duda, la más abundante y limpia. Esta es una razón importante, pero no la única, para

\* Corresponding author.

E-mail address: [mtausch@uni-wuppertal.de](mailto:mtausch@uni-wuppertal.de) (M.W. Tausch).

Peer Review under the responsibility of Universidad Nacional Autónoma de México.

Desplazamiento de Stokes;  
Conversión descendente;  
Efecto fotocromático;  
Solvatocromismo;  
Emisión inducida por agregación (AIE)

la inclusión de los fotoprocesos en la enseñanza de la química. En este artículo se presenta un enfoque experimental del concepto básico adecuado para la enseñanza de todos los fenómenos relacionados con la luz en una aproximación razonable. Se describen experimentos actualizados que muestran la conversión descendente de la luz por fluorescencia y fosforescencia a temperatura ambiente, así como la emisión inducida por la agregación y la luminiscencia inducida por la rigidificación. Con el fin de utilizar estos experimentos en la enseñanza de conceptos básicos en fotoquímica, se proponen una versión sencilla de un modelo de nivel de energía y una versión extendida, que contienen los estados vibracionales.

© 2017 Universidad Nacional Autónoma de México, Facultad de Química. Este es un artículo Open Access bajo la licencia CC BY-NC-ND (<http://creativecommons.org/licenses/by-nc-nd/4.0/>).

## Introduction

Light-involving phenomena are most suitable for the communication of core concepts in natural sciences in close combination with everyday experiences of students as well as with innovative technological applications. In this sense a series of teaching materials have been developed. Several of them are available *online* in English for free on our website (<http://www.chemiedidaktik.uni-wuppertal.de>; *Chemie und ihre Didaktik*, 2016a, 2016b) by clicking "Flash Animations" or "Teaching Photochemistry" respectively, further have been published in this Journal and others (Tausch, 2005; Tausch, Banerji, Scherf, 2013; Tausch & Bohrmann, 2003; Tausch, Bohrmann-Linde, Ibanez, et al., 2013; Tausch & Korn, 2001). Recent and good introductory texts regarding the importance of photochemistry, basic concepts and actual topics have been published by Bach (2015), Balzani, Bergamini, and Ceroni (2015), and Bléger and Hecht (2015).

According to N. J. Turro's paradigm of the "excited states of molecules" as "*the heart of all photoprocesses*" and the interpretation of the excited state as "*an electronic isomer of the ground state*" (Turro, 1978), it should be emphasized that the excited state  $A^*$  of a molecule A is not necessarily generated by light irradiation. It can also emerge from an exergonic chemical reaction or by supplying electrical energy (Fig. 1).


Actually, the excited state  $A^*$  can deactivate or react in very different ways as well. Therefore a veritable zoo of photoprocesses results from the different deactivation routes of  $A^*$  (Fig. 1). Traditionally some of them are allocated mainly to physics (i.e. photo- and electroluminescence). However, if one agrees with Turro that  $A^*$  is an electronic isomer of A, all photoprocesses have to be considered as photochemistry.

Due to the fact that fluorescence and phosphorescence phenomena nowadays count among everyday experiences of everyone, these should be used for teaching the core concepts of ground state and excited states in molecules.

## Experiments

### Phosphorescence at room temperature

Melt about 4 g of tartaric acid in a large test tube over the flame of a Bunsen burner. Heat gently to obtain a clear, molten fluid (it does not need to be colorless). Remove



**Figure 1** Generation and deactivation paths of the electronically excited state  $A^*$  (F: Fluorescence, P: Phosphorescence, CL: Chemiluminescence, EL: Electroluminescence, ECL: Electrochemiluminescence,  $S_1$ : 1<sup>st</sup> electronically excited singlet state,  $S_0$ : singlet ground state,  $T_1$ : 1<sup>st</sup> excited triplet state, D: Molecule of a compound D involved in energy transfer, B: Molecule of a compound B, involved in electron transfer, E: molecule of a starting compound E in a chemiluminescent reaction).

the heat and add about 5 mg of esculin to the molten substance. Shake slightly to dissolve the white powder in the tartaric acid and then spread the mixture on a large area of the test tube by tilting and rotating it (see video on <http://www.chemiedidaktik.uni-wuppertal.de/lehre/photo-mol/en/index.html>).

The fluorescence and phosphorescence effects can be explained by the model in Fig. 2, shown using common laboratory UV-light sources ( $\lambda = 366$  nm) or even a low cost (approx. 10\$) "UV-Flashlight" Ultrafire WF-501B,  $\lambda = 400$  nm, FWHM = 40 nm.

The fluorescence and phosphorescence spectra (Fig. 3) in this article have been recorded using a Varian Cary Eclipse Fluorescence Spectrometer by pouring the molten solution on a glass slide before rigidification.




Figure 2 The energy level model (basic version) – a core concept in chemistry.

### Temperature dependence of luminescence

Prepare three samples according to the instructions above. Heat one of the samples to 60°C using hot water. Cool another sample to 0°C using ice (or below in a freezer). Keep the third sample at room temperature. Irradiate all samples simultaneously and observe the color of the fluorescence and the occurrence and duration of the phosphorescence.

### Solvatochromism of spiropyran

Dissolve about 10 mg spiropiran in about 10 mL toluene, acetone and ethylene glycol respectively. Irradiate the solutions with a common laboratory UV-light source ( $\lambda = 366$  nm) or even a low cost "UV-Flashlight" (see above) for 60 s. Observe the different colors of the solutions and record spectra using a UV-VIS spectrometer.

The absorption spectra (Fig. 4) in this article have been recorded using a Analytik Jena AG Specord® 200 Plus UV VIS Spectrometer.



### Aggregation induced emission of merocyanine

Dissolve about 10 mg spiropiran in about 10 mL toluene and ethylene glycol respectively. Irradiate the solutions in a dark room with a common laboratory UV-light source ( $\lambda = 366$  nm) or even a low cost "UV-Flashlight" (see above). The bright red emission in ethylene glycol can be observed with the naked eye, while the toluene solution remains dark.

The corresponding fluorescence spectra (Fig. 4) in this article have been recorded using a Varian Cary Eclipse Fluorescence Spectrometer.

### Preparation of a rigid matrix of spiropiran in polystyrene

Dissolve about 30 mg spiropiran in 10 mL toluene (or xylene). Add about 3–4 g polystyrene until a viscous solution is obtained. The polystyrene should be styrofoam (i.e. chunks from packaging). Let the solution settle for a few minutes, so that the enclosed gas can leave the solution.



R = Glucose – Rest





Figure 3 Fluorescence and phosphorescence of esculin immobilized in tartric acid.



**Figure 4** Absorption and emission spectra of merocyanine in different solvents and in a solid matrix of polystyrene.

Fix a transparent foil to a level surface using tape. Spread the mixture thinly using a glass rod on the foil. Let the film dry for about 15 min in a fume cupboard and laminate it afterwards.

The rigidification induced fluorescence can then be shown using common laboratory UV-light sources ( $\lambda = 366 \text{ nm}$ ) or even a low cost "UV-Flashlight" (see above).

The corresponding isomerization is shown in Fig. 5, and the fluorescence spectrum is shown in Fig. 4. Fluorescence spectra in this article have been recorded using a Varian Cary Eclipse Fluorescence Spectrometer.

### Hazards and safety

Heeding the common laboratory safety regulation there are no hazards.


### Discussion and teaching recommendations

The color and the duration of light emission observable in the experiments described above provide convincing phenomena for introducing the concept of the ground state of molecules and the two principal types of excited states, the singlet state  $S_1$  and the triplet state  $T_1$ . In this sense the different colors of fluorescence and phosphorescence shown in Fig. 2 are in full agreement with theoretical predictions.


Similar experiments with fluorophors adsorbed on filter-paper (Schulman, 1976), immobilized in PMMA (Gilbert, 1977) or in boric acid (Roalstad, Rue, LeMaster, & Lasko, 1997) have been reported.

In contrast to some chemicals used in the references above, both tartric acid and esculin are non-hazardous and commercially available at a low price. Furthermore the preparation procedure of the samples as well as their investigation is fast and simple. We provide *online* on our website (<http://www.chemiedidaktik.uni-wuppertal.de/lehre/photo-mol/en/index.html>) a video clip showing the preparation of the luminescent sample and different model animations for explaining the Stokes shift of the fluorescence as well as the red shift of the phosphorescence compared to the fluorescence.

For teaching aims, at a *first* education level a very rough simplification of the Jablonski diagram is proposed in Fig. 2. The main reason for this simplification is the introduction of the concept of the ground state  $S_0$  and two types of excited states, the singlet state  $S_1$  and the triplet state  $T_1$ , both of them strongly related to observable and well distinguishable phenomena. In fact Fig. 2 provides theoretical reasons for the phenomenological differences between fluorescence and phosphorescence in our experiment, but it is not able to explain the Stokes shift – the red shift of fluorescence light compared to the exciting radiation. In our experiment however, the luminescent sample is irradiated with UV light at  $\lambda = 366 \text{ nm}$  and the maximum of the



**Figure 5** Photochromism: Spiropyran isomerizes reversibly into merocyanine.



**Figure 6** Energy level model (extended version) containing vibrational states and pathways of elementary processes in fluorescence and phosphorescence phenomena.

fluorescence emission occurs at  $\lambda = 421\text{ nm}$  as shown in Fig. 3. To explain this, an improved energy level model, including vibrational states inside of the same electronic states, is needed (Fig. 6). Without going into detail of the quantum mechanics involved (i.e. Born-Oppenheimer approximation and Franck-Condon principle), it can be made understandable that vibrational motions in molecules occur much more slowly (approx.  $10^{-12}\text{ s}$ ) than electron transitions (approx.  $10^{-15}\text{ s}$ ). In general the absorption of the exciting photon by a molecule leads to a higher vibrational level of the excited state  $S_1$  as shown in Fig. 6. The following vibrational relaxation (waved arrow in Fig. 6) occurs radiationless with the loss of a small amount of heat. From the lowest vibrational state of  $S_1$  the molecule deactivates by emitting a fluorescence photon, leading again to a higher vibrational state of  $S_0$ . From there it again relaxes into the lowest vibrational state of  $S_0$  by loss of heat. Following the complete cycle of elementary processes indicated in Fig. 6, one can understand why the emitted photon  $h\nu_2$  is less energetic than the absorbed photon  $h\nu_1$ .

This is the reason for the red shifted blue fluorescence compared to the exciting UV light (Stokes shift). In general fluorescence occurs as described here with a so called down conversion of photons.

Due to the fact that the lowest vibrational level of the excited triplet state  $T_1$  is even closer to the ground state, an additional down conversion of photons occurs by phosphorescence (Fig. 6). This is also a general rule and plausible given by the experimental observations and the used models described above.

According to the terminology proposed in literature (Mei et al., 2014), this is an example of *rigidification induced phosphorescence*. In the rigid matrix numerous intermolecular interactions, especially hydrogen bonds, collectively lock the conformations of the esculetin molecules. As a result of this the generation of excited triplet states and the deactivation via phosphorescence is improved.

The luminescence experiments at different temperatures show that the cooler the sample is, the longer the phosphorescence can be observed. The hot sample does not show

phosphorescence and the luminescence during irradiation changes from bluish-white to blue. These results confirm the transitions shown in Fig. 2, including the thermal equilibrium between  $S_1$  and  $T_1$ .


In order to demonstrate the phenomenon of *aggregation induced luminescence* (Ma et al., 2016; Mei et al., 2014) in classroom experiments, we use the hazard-free and commercially available spirobipyrrone, which is known for more than 60 years (Bergmann, Weizmann, & Fischer, 1950). Colorless spirobipyrrone in toluene (or xylene) solution is converted rapidly into its blue isomer merocyanine (Fig. 5) by irradiation with violet LED light or with sunlight (Tausch, 2005).

The reverse reaction of merocyanine to spirobipyrrone, observable by the decoloration of the solution, occurs spontaneously in darkness. It can be accelerated by heating or by irradiation with green light (Tausch, 2005; Tausch & Wachtendonk, 2007, <http://www.chemiedidaktik.uni-wuppertal.de/lehre/photo-mol/en/index.html>).

The results of the experiments regarding Aggregation Induced Emission (AIE) are summarized in Fig. 4.

As shown in Fig. 7, the color of merocyanine in toluene, acetone and ethylene glycol solution changes from blue over violet to red. Accordingly, the absorption maxima shift toward shorter wavelengths (negative solvatochromism; continuous lines in Fig. 4).

This is a characteristic behavior of dye molecules with a bipolar or even a zwitterionic resonance structure. Actually, merocyanine is an appropriate compound for explaining the phenomenon of negative solvatochromism. Again the energy level model is used, this time in an oversimplified version without vibrational levels (Fig. 7). The hypothesis is that, due to the electrostatic interaction of merocyanine zwitterions with the polar molecules of the solvent, the highest occupied energy level (HOMO) becomes lowered, whereas the lowest unoccupied energy level (LUMO) remains unmodified. Consequently, the energy  $\Delta E$  of the photons able to excite the merocyanine molecules in a more polar solvent has to be higher than in a less polar one. This effect should be increasingly distinct, the more polar the solvent



**Figure 7** Relation between the negative solvatochromism and the energies  $\Delta E$  needed for the electronic excitation of dye molecules in solvents of increasing polarities.

molecules are. Following this line of thought, one gets a pretty good idea for the behavior of merocyanine concerning light absorption in solvents of different polarities (Fig. 4, 7).

In the case of merocyanine the ability to emit light also depends strongly on the molecular environment. In our experiments, it is easily observable that a solution of merocyanine in the non-polar solvent toluene does not emit light, whereas a solution in the polar solvent ethylene glycol (as well as in a solid matrix of non-polar polystyrene) shows a bright red fluorescence.

The fluorescence of merocyanine in the polar solvent ethylene glycol and in the non-polar rigid matrix of polystyrene can and should be explained using the concepts of *aggregation induced luminescence* (Ma et al., 2016; Mei et al., 2014) and *rigidification induced luminescence* respectively (Yuan et al., 2010). In this sense one assumes that merocyanine zwitterion molecules aggregate in cages of polar solvent molecules. This aggregation leads to the restriction of intramolecular motions in merocyanine molecules. Thus the deactivation route of excited merocyanine molecules via vibronic relaxation is hindered or even suppressed and the deactivation by emission of photons is enhanced. In contrast to the ethylene glycol solution there is no aggregation of merocyanine molecules in the non-polar solvent toluene and intramolecular motions are not restricted. The deactivation of excited merocyanine molecules occurs mostly radiationless via vibronic relaxation. Immobilized in a rigid matrix, even in a non-polar one, the flexibility of merocyanine molecules becomes considerably limited and consequently the deactivation via fluorescence gets enhanced. Note that in both cases, in the polar ethylene glycol as well as in the non-polar polystyrene, according to the theoretical prediction (Fig. 6) the emitted photons are down converted compared to the absorbed photons. Actually, the absorption and emission spectra in Fig. 4 indicate a Stokes shift of approximately 130 nm.

## Conflicts of interest

The authors declare no conflicts of interest.

## Acknowledgment

We acknowledge the German Research Foundation (Deutsche Forschungsgemeinschaft, DFG) for supporting

the project "Photoprocesses in Science Education" ("Photoprozesse in der Lehre der Naturwissenschaften" Photo-LeNa, TA 228/4-1).

## References

- Bach, T. (2015). More chemistry with light! More light in chemistry!. *Angewandte Chemie International Edition*, 54(39), 11294–11295.
- Balzani, V., Bergamini, G., & Ceroni, P. (2015). Light: A very peculiar reactant and product. *Angewandte Chemie International Edition*, 54(39), 11320–11337.
- Bergmann, E. D., Weizmann, A., & Fischer, E. (1950). Structure and polarity of some polycyclic spirans. *Journal of the American Chemical Society*, 72(11), 5009–5012.
- Bléger, D., & Hecht, S. (2015). Visible-light-activated molecular switches. *Angewandte Chemie International Edition*, 54(39), 11338–11349.
- Chemie und ihre Didaktik, University of Wuppertal, <http://www.chemiedidaktik.uni-wuppertal.de/english/index.html> [consulted 20.11.16].
- Chemie und ihre Didaktik, University of Wuppertal: Photo-Mol, <http://www.chemiedidaktik.uni-wuppertal.de/lehre/photo-mol/en/index.html> [consulted 20.11.16].
- Gilbert, G. L. (1977). Phosphorescence: A demonstration. *Journal of Chemical Education*, 54(3), 167–168.
- Ma, X., Sun, R., Cheng, J., Liu, J., Gou, F., Xiang, H., et al. (2016). Fluorescence aggregation-caused quenching versus aggregation-induced emission: A visual teaching technology for undergraduate chemistry students. *Journal of Chemical Education*, 93(2), 345–350.
- Mei, J., Hong, Y., Lam, J. W. Y., Qin, A., Tang, Y., & Zhong Tang, B. (2014). Aggregation-induced emission: The whole is more brilliant than the parts. *Advanced Materials*, 26(31), 5429–5479.
- Roalstad, S., Rue, C., LeMaster, C. B., & Lasko, C. (1997). A room-temperature emission lifetime experiment for the physical chemistry laboratory. *Journal of Chemical Education*, 74(7), 853–854.
- Schulman, E. M. (1976). Room temperature phosphorescence – an experiment for the undergraduate physical or analytical laboratory. *Journal of Chemical Education*, 53(8), 522–524.
- Tausch, M. W. (2005). Teaching photochemistry with microscale experiments – a contribution for the ongoing curriculum modernization. *Educación Química*, 16(4), 529–533.
- Tausch, M. W., & Bohrmann, C. (2003). Photogalvanic cells for classroom investigations – a contribution for the ongoing curriculum modernization. *Journal of Chemical Education*, 80(12), 1471–1473.

- Tausch, M. W., Banerji, A., & Scherf, U. (2013). Classroom experiments and teaching materials on OLEDs with semiconducting polymers. *Educación Química*, 24(1), 17–22.
- Tausch, M. W., Bohrmann-Linde, C., Ibanez, J., Zavala-Araiza, D., Sotomayor-Martinez Barranco, B., Torres-Perez, J., et al. (2013). A demonstration of simultaneous electrochemoluminescence. *Journal of Chemical Education*, 90(4), 470–472.
- Tausch, M. W., & Korn, S. (2001). A laboratory simulation for coupled cycles of photosynthesis and respiration. *Journal of Chemical Education*, 78(9), 1238–1240.
- Tausch, M. W., & Wachtendonk, M. V. (Eds.). (2007). *CHEMIE 2000+, Textbook for High Schools*. Bamberg, Germany: C.C. Buchner (in German)
- Turro, N. J. (1978). *Modern molecular photochemistry*. New York, USA: Benjamin/Cummings.
- Yuan, W. Z., Shen, X. Y., Zhao, H., Lam, J. W. Y., Tang, L., Lu, P., et al. (2010). Crystallization-induced phosphorescence of pure organic luminogens at room temperature. *Journal of Physical Chemistry C*, 114(13), 6090–6099.