



Original article

Identifying molecularly defined antigens for a *Histoplasma capsulatum*-specific interferon gamma release assay



Marcela Rubio-Carrasquilla^{a,b}, Rodrigo Ochoa^c, Cristian Santa^d, Allan J. Guimarães^e, Luz Elena Cano^{a,f}, Ernesto Moreno^{g,*}

^a Grupo de Micología Médica y Experimental, Corporación para Investigaciones Biológicas (CIB), Medellín, Colombia

^b Instituto de Biología, Universidad de Antioquia, Medellín, Colombia

^c Programa de Estudio y Control de Enfermedades Tropicales – PECET, Universidad de Antioquia, Medellín, Colombia

^d Universidad Nacional de Colombia, Sede Medellín, Medellín, Colombia

^e Depto de Microbiología e Parasitología, Universidade Federal Fluminense, Niterói, RJ, Brazil

^f Escuela de Microbiología, Universidad de Antioquia, Medellín, Colombia

^g Facultad de Ciencias Básicas, Universidad de Medellín, Medellín, Colombia

ARTICLE INFO

Article history:

Received 30 November 2018

Accepted 14 June 2019

Available online 19 November 2019

Keywords:

Histoplasma capsulatum

IGRA

Molecularly defined antigens

Immunogenic proteins

T-cell epitopes

ABSTRACT

Background: In a previous work we showed the feasibility of an interferon gamma release assay (IGRA) for detecting latent infection by *Histoplasma capsulatum*. While in that proof-of-concept study we used crude fungal extracts as antigens, the newest IGRA developed for other infections are based on molecularly defined antigens, mostly on mixtures of immunogenic peptides.

Aims: To identify proteins in *H. capsulatum* that might serve as molecularly defined antigens for an IGRA test.

Methods: We surveyed the literature looking for known *H. capsulatum*-immunogenic proteins and assayed two of them as antigens in an IGRA test, in a study that involved 80 volunteers. Furthermore, we used several bioinformatics tools to identify specific *H. capsulatum* proteins and to analyze possible strategies for the design of *H. capsulatum*-specific immunogenic peptides.

Results: Seven *H. capsulatum*-immunogenic proteins were retrieved from the literature. IGRA tests using either the heat shock protein 60 or the M antigen showed high sensitivities but low specificities, most likely due to the high sequence similarity with the corresponding orthologs in other pathogenic microorganisms. We identified around 2000 *H. capsulatum*-specific proteins, most of which remain unannotated. Class II T-cell epitope predictions for a small number of these proteins showed a great variability among different alleles, prompting for a “brute force” approach for peptide design.

Conclusions: The *H. capsulatum* genome encodes a large number of distinctive proteins, which represent a valuable source of potential specific antigens for an IGRA test. Among them, the Cfp4 protein stands out as a very attractive candidate.

© 2019 Asociación Española de Micología. Published by Elsevier España, S.L.U. All rights reserved.

Identificación de antígenos definidos molecularmente para un ensayo de liberación de interferón-gamma específico para *Histoplasma capsulatum*

RESUMEN

Palabras clave:

Histoplasma capsulatum

IGRA

Antígenos definidos molecularmente

Proteínas inmunogénicas

Epítopes de células T

Antecedentes: En un trabajo anterior mostramos la viabilidad de un ensayo de liberación de interferón-gamma (IGRA) para detectar la infección latente por *Histoplasma capsulatum*. En esa prueba de concepto utilizamos extractos crudos del hongo como antígenos; sin embargo, los IGRA de última generación desarrollados para otras infecciones se basan en antígenos definidos molecularmente, principalmente en mezclas de péptidos inmunogénicos.

* Corresponding author.

E-mail address: emoreno@udem.edu.co (E. Moreno).

Objetivos: Identificar proteínas de *H. capsulatum* que podrían servir como antígenos definidos molecularmente en una prueba IGRA.

Métodos: Examinamos la literatura en busca de proteínas inmunogénicas de *H. capsulatum* ya conocidas, y ensayamos dos de ellas como antígenos en una prueba IGRA, en un estudio donde participaron 80 voluntarios. Además, utilizamos varias herramientas bioinformáticas para identificar proteínas específicas de *H. capsulatum* y analizar posibles estrategias para el diseño de péptidos inmunogénicos específicos.

Resultados: Encontramos siete proteínas de *H. capsulatum* caracterizadas como inmunogénicas en la literatura. Las pruebas IGRA donde utilizamos la proteína de choque térmico 60 o el antígeno M, mostraron una alta sensibilidad, pero baja especificidad, debido probablemente a la alta similitud de secuencia con los ortólogos correspondientes en otros microorganismos patógenos. Identificamos unas 2000 proteínas específicas de *H. capsulatum*, la mayoría de las cuales permanecen sin anotar. Las predicciones de epítopes de células T de clase II realizadas para un pequeño número de estas proteínas mostraron una gran variabilidad entre los diferentes alelos, sustentando la aplicación de un enfoque de «fuerza bruta» en el diseño de estos péptidos.

Conclusiones: El genoma de *H. capsulatum* codifica una gran cantidad de proteínas específicas que representan una fuente valiosa de posibles antígenos para una prueba IGRA. Entre ellos, la proteína Cfp4 resulta un candidato muy atractivo.

© 2019 Asociación Española de Micología. Publicado por Elsevier España, S.L.U. Todos los derechos reservados.

Histoplasmosis, caused by the dimorphic ascomycete fungus *Histoplasma capsulatum*, is the most common endemic mycosis in the Americas.⁹ Because immunocompromised individuals have a high risk of developing a disseminated *H. capsulatum* infection, the HIV pandemic and the use of immunosuppressive drugs have led to a dramatic increase in lethal histoplasmosis,^{2,23,28} developed either from primary infections and the re-burst of old ones with latent *H. capsulatum* foci.²⁵ Currently, no laboratory test can efficiently detect sub-clinical histoplasmosis, which may remain in a latent form even after several years from the primary infection.

In a previous proof-of-concept study²⁴ we showed the feasibility of an interferon gamma release assay (IGRA) specific for *H. capsulatum* to detect latent infection in asymptomatic individuals who suffered a *H. capsulatum* infection in the recent or distant past (from several months up to many years ago). The antigens tested in that first study were two histoplasmins and two yeast filtrates from different *H. capsulatum* strains (crude preparations).

In spite of the encouraging results obtained with the *H. capsulatum*-antigen crudes,²⁴ using these antigens may not be a suitable option for the development of a commercial IGRA test. Indeed, a survey of the literature shows that the different IGRA tests developed for several infection models have evolved from the use of crude antigens to the identification of unique proteins in the microorganism of interest and the use of peptide mixtures derived from these proteins. A relevant example of this is the evolution of QuantiFERON-TB for *Mycobacterium tuberculosis* (*Mtb*). Its first version was based on the same purified protein derivative (PPD) used in the traditional tuberculosis skin test (TST), whereas the next version (QuantiFERON-TB Gold) used a mixture of synthetic, overlapping peptides that fully cover the sequences of two *M. tuberculosis*-specific proteins (ESAT-6 and CFP-10), designed to be presented by class II HLA molecules and therefore being recognized by T-CD4+ lymphocytes.^{1,22}

Peptide mixtures following different conceptual designs have been used in several other IGRA tests targeting different infections, as cytomegalovirus (QuantiFERON-CMV),¹³ *Leishmania*,^{12,27} hepatitis B virus,⁸ varicella-zoster virus,²⁶ *Toxoplasma gondii*,⁷ *Mycobacterium leprae*²⁹ and HIV.¹⁴ Purified or recombinant proteins have been also used, as in *Leishmania* (histone B2 protein),²⁷ *M. leprae* (18, 65 and 70 kDa purified antigens)²⁹ and HIV-1 (p24 and gp120 recombinant proteins).¹⁶

In this work, using information from the literature, bioinformatics tools and assaying molecular candidates in an IGRA test that involved dozens of volunteers, we sought to identify proteins from *H. capsulatum* that might serve as molecularly defined antigens for a *H. capsulatum*-specific IGRA test, either in the form of whole recombinant antigens or as mixtures of synthetic peptides designed from their amino acid sequences.

Materials and methods

Identification of *H. capsulatum* antigens

Two alternative strategies were explored in parallel to identify possible molecularly defined antigens that might be used in a *H. capsulatum*-specific IGRA test: (i) Searching for *H. capsulatum* proteins that have been previously characterized as immunogenic molecules; and (ii) Identification of specific *H. capsulatum* proteins, that is, proteins that do not have close orthologs in other related fungi. In both cases, the identified proteins might be used in the form of complete proteins or as sources of specific immunogenic peptides. In the second approach, we used for comparison the sequences of all the proteins identified in the four *H. capsulatum* genomes that have been sequenced so far: G186AR, H88, H143 and NAm1, having 9315; 9449; 9582 and 18,736 identified proteins, respectively. Comparisons were made against the proteomes of four phylogenetically close fungal species: *Blastomyces dermatitidis* (11,211 proteins), *Emmonsia parva* (8562), *Paracoccidioides brasiliensis* (16,834) and *Paracoccidioides lutzii* (18,015). The identification of groups of orthologous proteins in these proteomes was carried out with the OrthoMCL program.²¹ Specific *H. capsulatum* proteins were analyzed with the TMHMM prediction method²⁰ to predict transmembrane domains; SignalP 3.0 and SecretomeP⁶ were used to predict signal peptides and secretion proteins, respectively, and TargetP¹¹ was used to predict mitochondrial proteins.

Sequences analyses and epitope prediction

The NCBI BLAST server (<https://blast.ncbi.nlm.nih.gov/Blast.cgi>) was used to compare the sequences of selected *H. capsulatum* proteins against the non-redundant NCBI sequence database. Predictions of class II epitopes were performed with the NetMHCII 2.3 server.¹⁹

Table 1

Sequence similarity (given as percentage of amino acid identity) among the seven immunogenic *H. capsulatum* proteins and their orthologs in related fungi and other unrelated microorganisms.

<i>H. capsulatum</i> immunogenic proteins	Genbank ID (gi)	Pathogenic microorganisms					
		<i>Emmonsia</i> <i>parva</i>	<i>Blastomyces</i> <i>dermatitidis</i>	<i>Paracoccidioides</i> <i>brasiliensis</i>	<i>Aspergillus</i> <i>fumigatus</i>	<i>Leishmania</i>	<i>Mycobacterium</i> <i>tuberculosis</i>
M antigen	225555166	75	76	62	65	— ^a	40 (67% cov.) ^b
HSP60	4558517	96	96	92	90	52	50
H antigen	225557355	86	87	80	67	—	34 (67% cov.)
H2B histone	32140475	96	96	95	94	43	—
HSP70	1230567	95	95	92	88	75	46
N-acetylated α-linked acidic dipeptidase	225557333	88	85	68	43	—	—
Cfp4	225555855	38 (51% cov.)	34 (57% cov.)	38 (52% cov.)	21 (63% cov.)	—	—

^a Sequence comparison did not yield significant results.

^b Percentage of the target sequence covered by BLAST. Not included if ≥85%.

Recombinant *H. capsulatum* antigens

The M antigen and the heat shock protein 60 (HSP60), both fused to a C-terminal His tag, were produced from previously obtained *Escherichia coli* clones.^{14,15} The proteins were purified by metal chelate chromatography using a Ni²⁺-Sepharose affinity column (TALON® Metal Affinity Resins; Clontech Laboratories Inc., USA). Protein yield was verified by SDS-PAGE, as well as by Western blot, using the serum from a *H. capsulatum*-infected individual as the source of primary antibodies against *H. capsulatum* antigens. Binding of serum antibodies was detected with an alkaline phosphatase-conjugated goat anti-human IgG antibody (Sigma-Aldrich, USA). For both antigens, thick bands corresponding to the expected molecular weights (between 72–95 kDa and ~60 kDa for the M antigen and HSP60, respectively) were observed.

IGRA tests

The study was conducted at the Corporación para Investigaciones Biológicas (CIB) in Medellin, Colombia. Eighty individuals who participated in our previous proof-of-concept study²⁴ volunteered for this second exploratory study. These were immunocompetent individuals with negative HIV test and age between 10 and 82 years. Pregnancy and immunosuppressive conditions were exclusive criteria. All the participants provided a written informed consent after receiving a thorough description of the study. Eighteen out of the 80 individuals never, presumably, had contact with *H. capsulatum* and were assigned to Group A (non *H. capsulatum*-infected), whereas 11 individuals (Group B-proven) showed clinical records proving that they were infected with *H. capsulatum* in the past. The remaining 51 individuals showed risk factors of exposure to *H. capsulatum* (Group B-risk). The assay was conducted as described by Rubio-Carrasquilla et al.²⁴ Briefly, whole blood (500 µl/well) was stimulated with 10 µg/ml of either HSP60 or M antigen. Non-stimulated blood-wells (basal level IFN-γ production) and blood-wells stimulated with 20 µg/ml of the T cell mitogen phytohemagglutinin (PHA) (Sigma-Aldrich, USA) were used as negative and positive controls, respectively. Cultures were incubated at 37 °C with 5% CO₂ for 24 h and then supernatants were collected. IFN-γ concentration in the supernatant was measured using a commercially available ELISA kit (Duo-set R&D Systems, Minneapolis, USA), following manufacturer's instructions. Optical density was measured with an ELISA microplate reader (Bio-Rad Laboratories Inc., USA) at 450 nm.

Statistical analyses

Statistical analyses were performed using R-project version 3.4.2. The Shapiro-Wilk normality test was employed to analyze the empirical distribution of the IFN-γ data, while the Levene's test was used to assess homoscedasticity. The Kruskal-Wallis non-parametric test and the post hoc Dunn's test were used for multiple comparisons of IFN-γ medians between the groups. A Benjamini-Hochberg correction was applied to control the type I error and *p*-adjusted (*q*) values were calculated. The assay sensitivity was evaluated as the percentage of positive responses in B-proven group, while specificity was calculated as the percentage of negative responses in the control group A. Statistical differences were considered significant if *p*-adjusted (*q*) was less than 0.05.

Ethics statement

All study procedures and written consent forms were approved by CIB's ethics review board.

Results and discussion

Immunogenic *H. capsulatum* proteins known from the literature

Table 1 lists seven proteins which have been evaluated before in multiple studies aiming to establish their role in the pathogenesis and immune response against *H. capsulatum*. Some of these proteins are part of the membrane and cell wall of the yeast phase, while others are secreted to the extracellular environment. Because of their proven immunogenicity, these proteins represent possible molecular antigens for a *H. capsulatum*-specific IGRA test.

IGRA tests using HSP60 or M antigen showed lack of specificity

We evaluated two of the *H. capsulatum* immunogenic proteins, HSP60 and M antigen, both in the form of recombinant proteins, as antigens in an IGRA test. The basal IFN-γ levels, in non-stimulated wells, were within a range of 0–36.8 pg/ml, with a median value of zero (data not shown). Fig. 1 shows the IFN-γ responses by group for each antigen. The two recombinant antigens were able to stimulate strong IFN-γ responses in the three groups. The response levels obtained for groups A (non *H. capsulatum*-infected) and B-proven show a strong overlapping for both antigens. For this reason, it is not possible to define a meaningful IFN-γ concentration cutoff value to discriminate between positive and negative responses. On the other




Fig. 1. Statistical box plots for IFN- γ production upon stimulation with PHA and the recombinant proteins HSP60 and antigen M. Symbols □, ■, ◆ represent individuals in groups A, B-proven and B-risk, respectively. The vertical axis is in logarithmic scale. For each individual, the IFN- γ response was adjusted by subtracting the corresponding basal level. IFN- γ concentrations between 0 and 0.1 pg/ml were assigned a value of 0.1 pg/ml to fit into the plot. Since the data have no normal distribution and are not homoscedastic, the Kruskal-Wallis non-parametric test and the post hoc Dunn's test were used for multiple comparisons of the medians (shown as thick bars within the boxes). The gray horizontal line is placed at the defined IFN- γ cutoff level (20 pg/ml).

hand we lack a gold standard to compare with, since none of the currently available tests is able to detect sub-clinical *H. capsulatum* infections. Therefore, we decided to use the cutoff value defined in our previous work, based on the IFN- γ responses to *H. capsulatum* crude antigens (20 pg/ml).²⁴

As expected from the strong overlapping of the IFN- γ responses, no statistically significant differences were found between groups A and B-proven for any of the antigens. Using the defined IFN- γ concentration cutoff, the obtained sensitivity values were high (91% and 81% for HSP60 and the M antigen, respectively), with very high IFN- γ response levels to the M antigen, in particular in the B-proven group. However, the specificities were low for both recombinant antigens (56% and 39%, respectively), with several members of group A (non *H. capsulatum*-infected) showing a medium-level or even a high IFN- γ response.

In B-risk group, 57% (29/51, for HSP60) or 63% (32/51, for M antigen) of the individuals were classified as infected using the above defined cutoff value. The percentage of positive responses in this group is much higher than the one obtained with the *H. capsulatum* crude antigens,²⁴ most likely because of the lower specificity shown by the two recombinant antigens. All the individuals in B-risk group had been classified as negative by the serology tests.²⁴ Thus, and in spite of its low specificity, the IGRA test performed here may be revealing a latent infection in a number of these individuals, especially in those yielding high (>100 pg/ml) IFN- γ responses.

It is worth recalling that M antigen (together with H antigen) is a standard protein tested by serology to detect *H. capsulatum* infection. Furthermore, both the M antigen and HSP60 are secreted by the fungus³ and therefore are components of the *H. capsulatum* crudes used as antigens in our first IGRA design, which showed a high specificity.²⁴ Possibly, using these antigens in pure form and, therefore, in greater quantities when compared to their fractions in the crude *H. capsulatum* preparations, should produce a higher level of antigen presentation, which might cause the exposure of epitopes that are cryptic in *H. capsulatum*, but dominant in other microorganisms.¹⁰ The results obtained indicate that neither HSP60 nor M antigen are suitable antigens for a *H. capsulatum*-specific IGRA test.

High sequence similarity between *H. capsulatum* immunogenic antigens and their orthologs in other pathogens

The sequences of the seven immunogenic *H. capsulatum* proteins were compared against those deposited in Genbank. Six of

these antigens showed a high similarity with their orthologs in related fungi (Table 1). The two heat shock proteins (HSP60 and HSP70), in addition, have similar orthologs in other less related pathogens such as *Leishmania* and *M. tuberculosis* (Table 1), as well as in other microorganisms (not shown). It is then plausible that the high IFN- γ responses produced by HSP60 and the M antigen in a significant number of non-*H. capsulatum* infected individuals are most likely cross reactions caused by the immune memory developed by these individuals against other infections. Furthermore, these results suggest that none of the first six proteins from Table 1 are suitable for an IGRA test, at least in the form of whole antigens.

Cfp4 is a *H. capsulatum*-specific protein

In contrast with the other six immunogenic proteins, Cfp4 has no close orthologs in related fungi. Furthermore, a BLAST search showed that there are no proteins similar to Cfp4 in any other organism (results not shown). In addition to its high specificity, Cfp4 is very attractive as a possible antigen for a *H. capsulatum*-specific IGRA test because of two other reasons: it is an immunogenic protein abundantly expressed and secreted in the yeast phase of the fungus¹⁷ and has a relatively small size, which allows covering its entire sequence with a limited number of overlapping peptides.

Prediction of *H. capsulatum*-specific proteins

We searched for proteins that are distinctive of *H. capsulatum* by comparing the sequences of all *H. capsulatum* proteins with the proteomes of other four related species of fungi, as detailed in Methods. The calculations performed with OrthoMCL yielded 11,344 homolog groups, out of which 1983 groups were composed only by homologous proteins from the four *H. capsulatum* strains. Since it is not practical to analyze each of these proteins individually, we focused on annotated proteins (that is, those with known or predicted function). This resulted in only 41 proteins. The next filter aimed to identify which of these proteins are secreted to the external medium. Using the SignalP server, we obtained five proteins that could be released by a classical secretion route. The remaining 36 proteins were processed with the SecretomeP server, which predicted 13 proteins associated with non-classical secretory pathways. The resulting 18 secretion proteins were further analyzed with the TargetP server, which indicated that none of them are mitochondrial proteins. These 18 proteins were also

Table 2

H. capsulatum-specific proteins selected from the orthology analysis.

Protein names	Genbank ID (gi)	Size (aa)
Aryl-alcohol dehydrogenase	225557852	105
Beta-lactamase	225558359	148
DUF636 domain-containing protein	225556157	272
Ankyrin repeat protein	225557463	268
Ankyrin repeat protein	225561435	529
Proteoglycan	225556461	666

analyzed with the TMHMM server, which predicted the presence of transmembrane helices in six of these proteins, that were consequently discarded.

BLAST server, the last filter, was used to compare the sequence of each of the remaining 12 proteins against the whole NCBI database of non-redundant proteins, to account for possible similarities with proteins from other pathogens, in addition to the four fungal species analyzed with OrthoMCL. Six out of the 12 proteins were discarded because of their significant sequence similarity (defined here as >50% amino acid identity in at least 70% coverage of the target sequence) with proteins from other pathogenic fungi, including species of *Emmonsia*, *Blastomyces*, *Sporothrix*, *Aspergillus*, *Fusarium*, *Trichophyton*, and *Coccidioides*, among others. Table 2 lists the six selected *H. capsulatum*-specific secretion proteins.

Cfp4 was not included within the 41 annotated proteins that were selected from the OrthoMCL output. The reason why Cfp4 is not annotated is precisely its lack of similarity with other proteins of known function. Cfp4 was, however, among the approximately 2000 proteins without close orthologs found by the program. Hence, there exists a vast space to be explored for possible specific antigens among the large number of non-annotated proteins found in the *H. capsulatum* genomes.

Peptide antigens from *H. capsulatum*-specific proteins: rational design or brute force?

As mentioned in the Introduction, the commercial IGAs available for *M. tuberculosis* and cytomegalovirus are based on peptide antigens. For *M. tuberculosis*, although the latest version (QuantiFERON-TB Plus) includes a second tube containing class I T-cell epitopes, the design of the test peptides is still aimed mainly at the stimulation of T-CD4+ lymphocytes.⁵ For *H. capsulatum*, a Th1-type response is needed for macrophage activation and efficient yeast clearance.¹⁸ In this type of response, T-CD4+ cells are crucial, being the main producers of IFN- γ .⁴ The T-CD8+ cells also play a role in host defense against histoplasmosis, but with minor consequences when compared to T-CD4+ lymphocytes, in

immunocompetent individuals.¹⁸ Therefore, a plausible approach toward a future development of a *H. capsulatum*-specific IGRA is the design of class II peptides from *H. capsulatum*-specific proteins, such as Cfp4 and the other distinctive proteins identified in this study.

We first explored a rational design approach of class II T-cell peptides. The NetMHCII-2.3 web server was employed to predict possible class II, 15 amino acid (aa) long epitopes in the seven distinctive *H. capsulatum* proteins identified here, for several HLA-DR, DP and DQ alleles serving as models for this exercise. Fig. 2 illustrates the results obtained for Cfp4 and beta-lactamase. Both proteins have a relatively small size, comparable to that of the *M. tuberculosis* proteins ESAT-6 and CFP-10. It can be observed that the number and distribution of predicted peptides throughout the sequences varies markedly from one allele to another. These results, together with the level of uncertainty that still characterizes the class II prediction algorithms, make the selection of a defined set of peptides (that should cover a wide spectrum of HLA alleles) risky. Taking also into account that the analyzed proteins are relatively short, so that their whole sequences can be covered with a limited number of overlapping peptides, a brute force approach seems more appropriate to create a pool of immunogenic peptides.

Conclusions

The development of an IGRA test for detecting latent infections by *H. capsulatum* should evolve toward the use of molecularly defined antigens, preferably immunogenic peptides. The use of known immunogenic *H. capsulatum* proteins having a high sequence similarity with their orthologs in other pathogenic microorganisms may result in a low test specificity, as demonstrated here for M and HSP60 antigens. Our bioinformatics analyses show that the *H. capsulatum* genome encodes a large number of distinctive proteins, which represent a valuable source of potential specific antigens. The Cfp4 protein, in particular, stands out as a very attractive candidate.

Funding

This work was supported by COLCIENCIAS (grant 221356933526) and the University of Antioquia (CODI-UdeA grant CIEMB-098-13).

Conflict of interest

The authors declare that they have no conflict of interest.




Fig. 2. Schematic representation of the epitopes predicted with the NetMHCII 2.3 server¹⁷ for different class II HLA alleles, for two distinctive proteins of *H. capsulatum*. The segments in gray represent 15 aa peptides with predicted affinities <500 nM. In cases where several consecutive peptides with the same nonapeptide nucleus were predicted, only the central peptide was represented.

References

1. Aagaard C, Brock I, Olsen A, Ottenhoff THM, Weldingh K, Andersen P. Mapping immune reactivity toward Rv2653 and Rv2654: two novel low-molecular-mass antigens found specifically in the *Mycobacterium tuberculosis* complex. *J Infect Dis.* 2004;189:812–9.
2. Adenis A, Aznar C, Couppié P. Histoplasmosis in HIV-infected patients: a review of new developments and remaining gaps. *Curr Trop Med Rep.* 2014;1:119–28.
3. Albuquerque PC, Nakayasu ES, Rodrigues ML, Frases S, Casadevall A, Zancopé-Oliveira RM, et al. Vesicular transport in *Histoplasma capsulatum*: an effective mechanism for trans-cell wall transfer of proteins and lipids in ascomycetes. *Cell Microbiol.* 2008;10:1695–710.
4. Allendoerfer R, Deepe GS. Intrapulmonary response to *Histoplasma capsulatum* in gamma interferon knockout mice. *Infect Immun.* 1997;65:2564–9.
5. Barcellini L, Borroni E, Brown J, Brunetti E, Codecasa L, Cugnata F, et al. First independent evaluation of QuantiFERON-TB Plus performance. *Eur Respir J.* 2016;47:1587–90.
6. Bendtsen JD, Jensen IJ, Blom N, Von Heijne G, Brunak S. Feature-based prediction of non-classical and leaderless protein secretion. *Protein Eng Des Sel.* 2004;17:349–56.
7. Chapey E, Wallon M, Debize G, Rabilloud M, Peyron F. Diagnosis of congenital toxoplasmosis by using a whole-blood gamma interferon release assay. *J Clin Microbiol.* 2010;48:41–5.
8. Dammermann W, Bentzien F, Stiel EM, Kühne C, Ullrich S, zur Wiesch JS, et al. Development of a novel IGRA assay to test T cell responsiveness to HBV antigens in whole blood of chronic hepatitis B patients. *J Transl Med.* 2015;13:157.
9. Deepe GS. *Histoplasma capsulatum* (histoplasmosis). In: Mandell, Douglas, and Bennett's principles and practice of infectious diseases; 2014. p. 2949–62.
10. Dominguez MR, Silveira EL, de Vasconcelos JR, de Alencar BC, Machado AV, Bruna-Romero O, et al. Subdominant/cryptic CD8 T cell epitopes contribute to resistance against experimental infection with a human protozoan parasite. *PLoS ONE.* 2011;6:1–12.
11. Emanuelsson O, Nielsen H, Brunak S, von Heijne G. Predicting subcellular localization of proteins based on their N-terminal amino acid sequence. *J Mol Biol.* 2000;300:1005–16.
12. Gidwani K, Jones S, Kumar R, Boelaert M, Sundar S. Interferon-Gamma release assay (modified quantiFERON) as a potential marker of infection for *Leishmania donovani*, a proof of concept study. *PLoS Negl Trop Dis.* 2011;5:e1042.
13. Giulieri S, Manuel O. QuantiFERON®-CMV assay for the assessment of cytomegalovirus cell-mediated immunity. *Expert Rev Mol Diagn.* 2011;11:17–25.
14. Guimaraes AJ, Hamilton AJ, de M Guedes HL, Nosanchuk JD, Zancopé-Oliveira RM. Biological function and molecular mapping of M antigen in yeast phase of *Histoplasma capsulatum*. *PLoS ONE.* 2008;3:e3449.
15. Guimaraes AJ, Frases S, Gomez FJ, Zancopé-Oliveira RM, Nosanchuk JD. Monoclonal antibodies to heat shock protein 60 alter the pathogenesis of *Histoplasma capsulatum*. *Infect Immun.* 2009;77:1357–67.
16. Hoher D, Benyoucef S, Chehadeh W, Hoher D, Benyoucef S, Chehadeh W, et al. Ex vivo interferon-gamma response to human immunodeficiency virus-1 derived peptides in human immunodeficiency virus-1 infected patients. *Scand J Immunol.* 2000;51:429–33.
17. Holbrook ED, Kemske MM, Richer SM, Wheat LJ, Rappleye CA. Glycosylation and immunoreactivity of the *Histoplasma capsulatum* Cfp4 yeast-phase exoantigen. *Infect Immun.* 2014;82:4414–25.
18. Horwath MC, Fecher RA, Deepe GS. *Histoplasma capsulatum*, lung infection and immunity. *Future Microbiol.* 2015;10:967–75.
19. Jensen KK, Andreatta M, Marcattili P, Buus S, Greenbaum JA, Yan Z, et al. Improved methods for predicting peptide binding affinity to MHC class II molecules. *Immunology.* 2018;154:394–406.
20. Krogh A, Larsson B, Von Heijne G, Sonnhammer ELL. Predicting transmembrane protein topology with a hidden Markov model: application to complete genomes. *J Mol Biol.* 2001;305:567–80.
21. Li L, Stoeckert CJ, Roos DS. OrthoMCL: identification of ortholog groups for eukaryotic genomes. *Genome Res.* 2003;13:2178–89.
22. Masurek GH, Villarino ME. Guidelines for using the QuantiFERON-TB test for diagnosing latent *Mycobacterium tuberculosis* infection. Centers for Disease Control and Prevention. *MMWR Recomm Rep.* 2003;52:15–8.
23. Nel JS, Bartelt LA, van Duin D, Lachiewicz AM. Endemic mycoses in solid organ transplant recipients. *Infect Dis Clin N Am.* 2018;32:667–85.
24. Rubio-Carrasquilla M, Santa CD, Rendón JP, Botero-Garcés J, Moreno E, Cano LE. An interferon gamma release assay specific for *Histoplasma capsulatum* to detect asymptomatic infected individuals: a proof of concept study. *Med Mycol.* 2019;57:724–32.
25. Sifuentes-Osornio J, Corzo-León DE, Ponce-De-León IA. Epidemiology of invasive fungal infections in Latin America. *Curr Fungal Infect Rep.* 2012;6:23–34.
26. Smith JG, Liu X, Kaufhold RM, Clair J, Caulfield MJ. Development and validation of a gamma interferon ELISPOT assay for quantitation of cellular immune responses to varicella-zoster virus. *Clin Diagn Lab Immunol.* 2001;8:871–9.
27. Turgay N, Balcioğlu IC, Toz SO, Ozbel Y, Jones SL. Short report: Quantiferon-*Leishmania* as an epidemiological tool for evaluating the exposure to *Leishmania* infection. *Am J Trop Med Hyg.* 2010;83:822–4.
28. Vergidis P, Avery RK, Wheat LJ, Dotson JL, Assi MA, Antoun SA, et al. Histoplasmosis complicating tumor necrosis factor- α blocker therapy: a retrospective analysis of 98 cases. *Clin Infect Dis.* 2015;61:409–17.
29. Weir RE, Morgan AR, Britton WJ, Butlin CR, Dockrell HM. Development of a whole blood assay to measure T cell responses to leprosy: a new tool for immuno-epidemiological field studies of leprosy immunity. *J Immunol Methods.* 1994;176:93–101.